

KellyMed

Blood and Infusion Warmer



5"
Touch Screen

Double
Channel

Microprocessor
Control

KL-2031N



Innovative Infusion Systems, Since 1994

Beijing KellyMed Co., Ltd.

Blood and Infusion Warmer

KL-2031N

Technical Specifications

Product name	Blood and Infusion Warmer
Model	KL-2031N
Application	Blood transfusion, infusion, enteral nutrition, parenteral nutrition,
Warmer Channel	Double Channel
LCD	5" Touch screen
Temperature	30-42°C, in 0.1°C increments
Temperature accuracy	±0.5°C
Warm time	<3 mins from 23±2°C to 36°C
Alarms	Over temperature alarm, low temperature alarm warm malfunction, low battery
Additional Features	Real-time temperature, automatic power switching, programmable fluid name and temperature range.
Wireless Management	Optional
Power Supply, AC	100-240 V, 50/60 Hz, ≤100 VA
Battery	18.5 V, rechargeable
Battery Life	5 hours for single channel, 2.5 hours for double channel
Working Temperature	0-40°C
Relative Humidity	10-90%
Atmospheric Pressure	860-1060 hpa
Size	110(L)*50(W)*195(H) mm
Weight	0.67 kg
Safety Classification	Class II, type CF
Fluid Ingress Protection	IP43

Beijing KellyMed Co., Ltd.

Add: 6R International Metro Center, No. 3 Shilipu,
Chaoyang District, Beijing, 100025, China

Mobile/WhatsApp/WeChat: +86-13671114111

E-mail: export@kelly-med.com

Web: www.kelly-med.com