

Precision filtration infusion set for pumps


Our Technology Your Health



- Special for infusion pumps
- Precision filtration
- Multifunctional care
- Efficient filtration
- DEHP free

Jiangsu JevKev Medtec Co.,Ltd.

Application scope



Regulatory requirements

- The national standards, registration certificates, user manuals, and packaging of gravity infusion sets clearly require that they are only for gravity use and should not be used for pressure infusion.
- To distinguish between pump infusion and gravity infusion, use a matching pump infusion set when using an infusion pump.
- When using infusion pumps internationally, they are equipped with infusion sets.

Cork puncture device

Thick handle design, labor-saving, sturdy and resilient.



Standard pore size of filter medium

2um、3um、5um



Regulator



Roller slide switch, convenient control.

Functional advantages

- Filter aperture: 2um、3um、5um.
- Characteristics: Super hydrophilic, good mechanical properties, no fiber shedding, suitable for any environment.
- Double membrane: exhaust membrane+ antibacterial membrane.
- Material: Special medicinal liquid filtration material, with a filtration rate of 99%, effectively reducing the occurrence of infusion reactions.


Special for infusion pumps

- The inner cavity wall is smooth, non adhesive after prolonged compression, resistant to fatigue, high tear resistance, and high rebound.
- Pressure resistant tubing, capable of withstanding a pressure of 200KPA.
- Non phthalic material, DEHP free, protects patients and medical staff from Plasticizer hazards.



Security

- It can effectively block particles from entering the human body.
- Reduce the occurrence of adverse reactions such as phlebitis, local ischemic edema, and thrombosis.
- Prevent the occurrence of heat source like reactions.



KellyMed | **Jev&Kev**



Jiangsu JevKev Medtec Co.,Ltd.

Address: East 1st & 2nd Floor, Building G06, No.1 Avenue China Medical City, TaiZhou Medical Hi-Tech Industrial Development Zone, 225300, Taizhou city, Jiangsu Province, China.

Tel: +86-523-86050096

FAX: +86-523-8681 0056

WEB: www.jevkev.com

E-MAIL: international@kelly-med.com

SINCE 1994

Leading the safe infusion for 30 years